



Recombinant Human CREB3 regulatory factor (CREBRF)

Product Code	CSB-MP815547HU
Abbreviation	CREBRF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUR6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPQPSVSGMD PPFGDAFRSH TFSEQTLMST DLLANSSDPD FMYELDREMNYQQNPRDNFL SLEDCKDIEN LESFTDVLDN EGALTSNWEQ WDTYCEDLTK YTKLTSCDIW GTKEVDYLGL DDFSSPYQDE EVISKTPTLA QLNSEDSQSV SDSLYYPDSL FSVKQNPLPS SFPGKKITSR AAPVCSSKT LQAEVPLSDC VQKASKPTSS TQIMVKTNMY HNEKVNHFVE CKDYVKKAKV KINPVQQSRP LLSQIHTDAA KENTCYCGAV AKRQEKKGME PLQGHATPAL PFKETQELL SPLPQEGPGS LAAGESSSLASTSVSDSSQ KKEEHNYSLF VSDNLGEQPT KCSPEEDEED EEDVDDDEDHD EGFGSEHELSENEEEEEEEEE DYEDDKDDDI SDTFSEPGYE NDSVEDLKEV TSISSRKRKG RRYFWEYSEQ LTPSQQERML RPSEWNRDTL PSNMYQKNGL HHGKYAVKKS RRTDVEDLTP NPKKLLQIGN ELRKLNKVIS DLTPVSELPL TARPRSRKEK NKLASRACRL KKKAQYEANK VKLWGLNTEY DNLLFVINSI KQEIVNRVQN PRDERGPNMG QKLEILIKDT LGLPVAGQTS EFVNQVLEKT AEGNPTGGLV GLRIPTSKV
Source	Mammalian cell
Target Names	CREBRF
Protein Names	Recommended name: CREB3 regulatory factor Alternative name(s): Luman recruitment factor Short name= LRF
Expression Region	1-639
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.