



Recombinant Human Cyclin-J-like protein (CCNJL)

Product Code	CSB-MP815550HU
Abbreviation	CCNJL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IV13
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMDEPWWEGR VASDVHCTLR EKELKLPTFR AHSPLLKSRR FFVDILTLLS SHCQLCPAAR HLAVYLLDHF MDRYNVTTSK QLYTVAVSCL LLANGVSLLS PRLKCSGMIS AHCNLHLP GS SNPASAPHP PPTPPQVAET TGKFEDREDH VPKLEQINST RILSSQNFTL TKKELLSTEL LLEAFSWNL CLPTPAHFLD YYLLASVSQK DHHCHTWPTT CPRKTKECLK EYAHYFLEV T LQDHIFYKFQ PSVVAACV G ASRICLQLSP YWTRDLQRIS SYSLEHLSTC IEILLVVYDN VLKDAVAVKS QALAMVPGTP PTPTQVLFQP PAYPALGQPA T TLAQFQTPV QDLCLAYRDS LQAHRSGSLL SGSTGSSLHT PYQPLQPLDM CPVVPASLS MHMAIAAEPR HCLATTYGSS YFSGSHMFPT GCFDR
Source	Mammalian cell
Target Names	CCNJL
Protein Names	Recommended name: Cyclin-J-like protein
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.