



Recombinant Human Ankyrin repeat and MYND domain-containing protein 2 (ANKMY2)

Product Code	CSB-MP815551HU
Abbreviation	ANKMY2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IV38
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVHIKKGELT QEEKELLEVI GKGTVQEAGT LLSSKNVRVN CLDENGMTPL MHAAYK GKLD MCKLLLRHGA DVNCHQHEHG YTALMFAALS GNKDITWVML EAGAETDVVN SVGR TAAQMA AFVGQHDCVT IINFFPRER LDYYTKPQGL DKEPKLPPKL AGPLHKIITT TNLHPVKIVM LVNENPLLTE EAALNKCYRV MDLICEKCMK QRDMNEVLAM KMHYISCIFQ KCINFLKDGE NKLDTLIKSL LKGRASDGFP VYQEKIRES IRKFPYCEAT LLQQLVRSIA PVEIGSDPTA FSVLTQAITG QVGFVDVEFC TTCGEKGASK RCSVCKMVIY CDQTCQKTHW FTHKKICKNL KDIYEKQQLE AAKEKRQEEN HGKLDVNSNC VNEEQPEAEV GISQKDSNPE DSGEGKKESL ESEAELEGLQ DAPAGPQVSE E
Source	Mammalian cell
Target Names	ANKMY2
Protein Names	Recommended name: Ankyrin repeat and MYND domain-containing protein 2
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.