



# Recombinant Human C-Maf-inducing protein (CMIP)

<b>Product Code</b>	CSB-EP815577HU
<b>Abbreviation</b>	CMIP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IY22
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDVTSSSSGGG GDPRQIEETK PLLGGDVSAP EGTKMGAVPC RRALLCNGM RYKLLQEGDI QVCVIRHPRT FLSKILTSKF LRRWEPHHLT LADNSLASAT PTGYMENSVS YSAIEDVQLL SWENAPKYCL QLTIPGGTVL LQAANSYLRD QWFHSLQWKK KIYKYKKVLS NPSRWEVVLK EIRTLVDMAL TSPLQDD SIN QAPLEIVSKL LSENTNLTTQ EHENIIVAIA PLENNHPPP DLCEFFCKHC RERPRSMVVI EVFTPVVQRI LKHNMDFGKC PRLRLFTQEY ILALNELNAG MEVVKKFIQS MHGPTGHCPH PRVLPNLVAV CLAAIYSCYE EFINSRDNSP SLKEIRNGCQ QPCDRKPTLP LRLHPSPDL VSQEATLSEA RLKSVVVASS EIHVEVERTS TAKPALTASA GNDSEPNLID CLMVSPACST MSIELGPQAD RTLGCYVEIL KLLSDYDDWR PSLASLLQPI PFPKEALAHE KFTKELKYVI QRFAEDPRQE VHSCLLSVRA GKDGWFQLYS PGGVACDDDG ELFASMVHIL MGSCYKTKKF LLSLAENKLG PCMLLALRGN QTMVEILCLM LEYNIIDNND TQLQIISTLE STDVGKRMYE QLCDRQRELK ELQRKGGPTR LTLPSKSTDA DLARLLSSGS FGNLENLSLA FTNVTSAEAE HLIKLP SLKQ LNLWSTQFGD AGLRLLSEHL TMLQVLNLCE TPVTDAGLLA LSSMKSLCSL NMNSTKLSAD TYEDLKAKLP NLKEVDVRYT EAW
<b>Source</b>	E.coli
<b>Target Names</b>	CMIP
<b>Protein Names</b>	Recommended name: C-Maf-inducing protein Short name= c-Mip Alternative name(s): Truncated c-Maf-inducing protein Short name= Tc-Mip
<b>Expression Region</b>	1-773
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.