



Recombinant Human C-Maf-inducing protein (CMIP)

Product Code	CSB-EP815577HU
Abbreviation	CMIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IY22
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDVTSSSGGG GDPRQIEETK PLLGGDVSAP EGTKMGAVPC RRALLCNGM RYKLLQEGDI QVCVIRHPRT FLSKILTSKF LRRWEPHHLT LADNSLASAT PTGYMENSVS YSAIEDVQLL SWENAPKYCL QLTIPGGTVL LQAANSYLRD QWFHSLQWKK KIYKYKKVLS NPSRWEVVLK EIRTLVDMAL TSPLQDD SIN QAPLEIVSKL LSENTNLTTQ EHENIIVAIA PLENNHPPP DLCEFFCKHC RERPRSMVVI EVFTPVVQRI LKHNMDFGKC PRLRLFTQEY ILALNELNAG MEVVKKFIQS MHGPTGHCPH PRVLPNLVAV CLAAIYSCYE EFINSRDNSP SLKEIRNGCQ QPCDRKPTLP LRLHPSPDL VSQEATLSEA RLKSVV VASS EIHVEVERTS TAKPALTASA GNDSEPNLID CLMVSPACST MSIELGPQAD RTLGCYVEIL KLLSDYDDWR PSLASLLQPI PFPKEALAHE KFTKELKYVI QRFAEDPRQE VHSCLLSVRA GKDGWFQLYS PGGVACDDDG ELFASMVHIL MGSCYKTKKF LLSLAENKLG PCMLLALRGN QTMVEILCLM LEYNIIDNND TQLQIISTLE STDVGKRMYE QLCDRQRELK ELQRKGGPTR LTLPSKSTDA DLARLLSSGS FGNLENLSLA FTNVTSACAE HLIKLP SLKQ LNLWSTQFGD AGLRLLSEHL TMLQVLNLCE TPVTDAGLLA LSSMKS L CSL NMNSTKLSAD TYEDLKAKLP NLKEVDVRYT EAW
Source	E.coli
Target Names	CMIP
Protein Names	Recommended name: C-Maf-inducing protein Short name= c-Mip Alternative name(s): Truncated c-Maf-inducing protein Short name= Tc-Mip
Expression Region	1-773
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.