



# Recombinant Human Coiled-coil domain-containing protein C6orf97 (C6orf97)

<b>Product Code</b>	CSB-BP815585HU
<b>Abbreviation</b>	C6orf97
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYT3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLDCTSHIA LGAASPAPPEE TYDHLSEVPV TREQLNHYRN VAQNARSELA ATLVKFECAQ SELQDLRSKM LSKEVSCQEL KAEMESYKEN NARKSSLLTS LRDRVQELEE ESAALSTSKI RTEITAAHAI KENQELKKKV VELNEKLQKC SKENEENKKQ VSKNCRKHEE FLTQLRDCLD PDERNDKASD EDLILKLRDL RKENEFVKGQ IVILEETINV HEMEAkasRE TIMRLASEVN REQKKAASCT EEKEKLNQDL LSAVEAKEAL EREVKIFQER LLAGQQVWDA SKQEVSLKK SSSELEKSLK ASQDAVTTSQ SQYFSFREKI AALLRGRISM TGSTEDTILE KIREMDSREE SRDRMVSQLE AQISELVEQL GKESGFHQKA LQRAQKAENM LETLQGQLTH LEAELVSGGV LRDNLNFEKQ KYLKFLDQLS QKMKLDQMAA ELGFDMRLDV VLARTEQLVR LESNAVIENK TIAHNLQRKL KTQKERLESK ELHMSLLRQK IAQLEEEKQA RTALVVERDN AHLTIRNLQK KVERLQKELN TCRDLHTELK AKLADTNELK IKTLEQTKAI EDLNKSRDQL EKMKEKAEKK LMSVKSELDT TEHEAKENKE RARNMIEVVT SEMKTLKKSLE EAEKREKQL ADFREVVSQM LGLNVTSLAL PDYEIHKCLE RLVHSHQHFF VTCACLKDV TGQERHPQGH LQLLH
<b>Source</b>	Baculovirus
<b>Target Names</b>	CCDC170
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein C6orf97
<b>Expression Region</b>	1-715
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.