



Recombinant Human Protein FAM9A (FAM9A)

Product Code	CSB-EP815595HU-B
Abbreviation	FAM9A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZU1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPVGRKRSR KAAKAQLEAQ VTA AQGATKE GSGIASNFPG QPTMEPVGRK RSRKA AKAQL EAQVRAAPAK KHTGKDPVRD ECEERNPFTE TREEDVTDEH GEREPFAEKD EHTGIHTMKL EHIAADIKKG LAAKREMIKI DKAAYRKTKN TIERALKKKQ LKRQKR DYRH TRKLLNVLKE YIAEKQKDDE AEEAEAAAAA AEAAAAAEAA AAAAEVIVVE DEEEEEKEEEE EEKEEEEEEG EEEGGGEEGE EGGGGGEGEE TEEEEEEEEEE EEEEEQIKAF QEKQKRWQQP TGVR SWRLRE MKPLLEQLLK AAKDTKDNYC IISSESESEL DN
Source	E.coli
Target Names	FAM9A
Protein Names	Recommended name: Protein FAM9A
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.