



Recombinant Mouse RNA-binding protein 41 (Rbm41)

Product Code	CSB-EP815630MO
Abbreviation	Rbm41
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8JZV4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MKRVNSCVKD EEHVLEELET EGERQLKSLL QHQLDTSVSI EECVSKKKSF APGTMYPFG KEAAGTMTLS QFQTLHEKDQ ETASLRELGL NETEILWKS HVSGEKRTKL RATPEAIQKR LEDIKERISE RQRILCLPQR FSKSKQLTRR EMEIEKSLFQ GTDRHSFLKA LYYQDEPPKK NKGDPMNNLE HFYRETIMKK RLEEFQLLRG ESFACHSLVS AASVSGSGTA EKPSLLQDKG KQAAQGGKPR LHVAKLIDFP TEQYWTGPKT LKQPIEFIFE DEIQNRNLSE EEIRNIPMFS SYNPGEPNKV LYLKLNLSRV KERDLISLFA RFQEKKGPP I QFRMMTGRMR GQAFLTFPNK DIAWQALHQI NGYKLYGKIL VIEFAKSKKQ QSAVQRSSLI TSAPDGRTEI NGS
Source	E.coli
Target Names	Rbm41
Protein Names	Recommended name: RNA-binding protein 41 Alternative name(s): RNA-binding motif protein 41
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.