



# Recombinant Mouse Protein OS-9 (Os9)

<b>Product Code</b>	CSB-BP815653MO
<b>Abbreviation</b>	Os9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K2C7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GVGS LNLEELSEMR YGIQILPLPV MGGQSQASDV VVVSSKYKQR YECRLPAGAI HFQREREET PAYQGPGIPE LLSPMRDAPC LLKTKDWWTY EFCYGRHIQQ YHMEDSEIKG DVLYLGHYQS SFNWDEETAK ASKQHRLKRY HSQTYGNGSK CDLNGKPREA EVRFLCDEGA GISGDYIDRV DEPVSCSYVL TIRTSRLCPH PLLRPPASAA PQAILCHPAL QPDEYMAYLQ RQAESKQHEE KTTEEVQDTD RQVWSGSKAA GAPPKKEDVS PAKEEKESEL WKLQGPEEQA AAREEAQAGE QDLNHEAAD PAPSPPNDFQ NNVQVKLIRS PADLIRLIEE LKAAEKGKPS VRREQPGDDT TEAPQREAEG TKAKGKDGEF PGLMEEEDGD DEEEEEEEEE DEEEQQLLGE FEKELEGMLL PSNRERLRSE VKAGMERELE NIIQETEKEL DPEGLRKESE REQAILALTS TLDKLIKRLQ ENQSPELVQK YKKRRVVPQK PPPSPHPTTE EPEHRVRVRV TKLRPGGPNR DLTVLEMNRE NPQLKQIEGL VTEVLEREGL TAEGKIEIKI VRPGAEGKEE DTRWLTDEDT RNLKEIFFNI LVQGAEANK ERQRQSELES NYRRVWGSPG GEDTGDLDEF DF
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os9
<b>Protein Names</b>	Recommended name: Protein OS-9
<b>Expression Region</b>	27-672
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.