



Recombinant Mouse Protein OS-9 (Os9)

Product Code	CSB-BP815653MO
Abbreviation	Os9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K2C7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	GVGS LNLEELSEMR YGIQILPLPV MGGQSQASDV VVVSSKYKQR YECRLPAGAI HFQREREET PAYQGPGIPE LLSPMRDAPC LLKTKDWWTY EFCYGRHIQQ YHMEDSEIKG DVLYLGHYQS SFNWDEETAK ASKQHRLKRY HSQTYGNGSK CDLNGKPREA EVRFLCDEGA GISGDYIDRV DEPVSCSYVL TIRTSRLCPH PLLRPPASAA PQAILCHPAL QPDEYMayLQ RQAESKQHEE KTTEEVQDTD RQVWSGSKAA GAPPKKEDVS PAKEEKESEL WKLQGPEEQA AAREEAQAGE QDLNHEAAAD PAPSPPNDFQ NNVQVKLIRS PADLIRLIEE LKAAEKGKPS VRREQPGDDT TEAPQREAEG TKAKGKDGEF PGLMEEEDGD DEEEEEEEEE DEEEQQLLGE FEKELEGMLL PSNRERLRSE VKAGMERELE NIIQETEKEL DPEGLRKESE REQAILALTS TLDKLIKRLQ ENQSPELVQK YKKRRVVPQK PPPSPHPTEE EPEHRVRVRV TKLRPGGPNR DLTVLEMNRE NPQLKQIEGL VTEVLEREGL TAEGKIEIKI VRPGAEGKEE DTRWLTDEDT RNLKEIFFNI LVQGAEANK ERQRQSELES NYRRVWGSPG GEDTGDLEDF DF
Source	Baculovirus
Target Names	Os9
Protein Names	Recommended name: Protein OS-9
Expression Region	27-672
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.