



Recombinant *Oryza sativa* subsp. *japonica* DEAD-box ATP-dependent RNA helicase 36 (Os07g0633500, LOC_Os07g43980), partial

Product Code	CSB-EP815810OFG-B
Abbreviation	Os07g0633500, LOC_Os07g43980
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L4E9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MEVDGEARPF LLFSKPKSSK KKPQKQEAEPQ VHTQPEEPPN PSPSPAIEPD LRDSDEAPAA AVTEHAGDDA AAAAVPSTFA ELGLSQWLVD VCDLGMVRVP TAVQRRICIPR ALEGRDVLGI AETGSGKTAA FALPILHRLG EDPYGVAALA LAPTRELAQ LAEQFRALGA PLGLRCLAAI GGFDSLQAK GLARRPHVVV ATPGRIATLI NDDPDLAKVF ARTKFLVLDE ADRVLDINFE EDLRVIFGSL PKKRQTFIFS ATISDNLRSL LELSGNNSYF FEAYEGFKTV DTLKQLYIHV PPDAKELYLF YLLSKMNEDN IRSVIVFVST CRTCYLDLFL LEELGHPAVS LHSHKPQSRR LAALHNFKSS KVPVLLATDV ASRGLDIQTV DLVINYDVPR YPRDYIHRVG RTARATRGGL SISFITTQRD IRLLEIEDV VGKQLGAYDG EMRDVNKDAT KVFKAARLAN MKMADEGHED KVQARKEQKK RAQERKRKHD E
Source	E.coli
Target Names	Os07g0633500
Protein Names	Recommended name: DEAD-box ATP-dependent RNA helicase 36 EC=3.6.4.13
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.