



Recombinant Arabidopsis thaliana WD repeat-containing protein DWA1 (DWA1)

Product Code	CSB-BP815811DOA
Abbreviation	DWA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L4M1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MYGDATNWE DEYRESILKE REIETRTVFR TAWAPPARIS NPDAFVVASS DGTlafHSLN SLVSQSASFG YSKGQDVMVA EPERVVRAHE GPAYDVKFGY EDEDALLSC GDDGRVRGWK WREFAESDVS LHLKENHLKP LLELINPQHK GPWGALSPMP EINAMSVDPQ SGSVFtaAGD SCAYCWDVES GKIKMTFKGH SDYLHTVVSr SSASQILTGS EDGTARIWDC KTGKCVKVIg SQDKKSRLRV SSMALDGSES WLVCgQGKNL ALWNLpASEC VQTIPIAHV QDVMFDEKQI LTVGAEPllR RFDLNgALLS QIHcAPCSVF SISLHPAGVV AVGGYGGIVD VISQFGSHLC TFRSSSL
Source	Baculovirus
Target Names	THO6
Protein Names	Recommended name: WD repeat-containing protein DWA1 Alternative name(s): Protein DWD HYPERSENSITIVE TO ABA 1
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.