



Recombinant Mouse Serine/threonine-protein kinase PAK 2 (Pak2)

Product Code	CSB-MP815882MO
Abbreviation	Pak2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CIN4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SDNGELEDK PPAPPVRMSS TIFSTGGKDP LSAHNSLKPL PSVPPEKKPR NKIISIFSGT EKGSKKKEKE RPEISPPSDF EHTIHVGFDA VTGEFTGMPE QWARLLQTSN ITKLEQKKNP QAVLDVLKFY DSNTVKQKYL SFTPPEKDG PSGTPALNTK GSETSAVVTE EDDDDEDAAP PVIAPRPDHT KSIYTRSVI PIPAPVGDSN VD SGAKSSDK QKKKAKMTDE EIMEKLRTIV SIGDPKKKYT RYEKIGQGAS GTVFTATDVA LGQEVAIKQI NLQKQPKKEL IINEILVMKE LKNPNIVNFL DSYLVGDEL VVMEYLAGGS LTDVVTETCM DEAQIAAVCR ECLQALEFLH ANQVIHRDIK SDNVLGMEG SVKLDFGFC AQITPEQSKR STMVGTPYWM APEVVTRKAY GPKVDIWSLG IMAIEMVEGE PPYLNENPLR ALYLIATNGT PELQNPEKLS PIFRDFLNRC LEMDVEKRG S AKELLQHPFL KLAKPLSSLT PLILAAKEAM KSNR
Source	Mammalian cell
Target Names	Pak2
Protein Names	Recommended name: Serine/threonine-protein kinase PAK 2 EC= 2.7.11.1 Alternative name(s): Gamma-PAK p21-activated kinase 2 Short name= PAK-2 Cleaved into the following 2 chains: 1. PAK-2p27 2. PAK-2p34
Expression Region	2-524
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The p21 activated kinases (PAK) are critical effectors that link Rho GTPases to cytoskeleton reorganization and nuclear signaling. The PAK proteins are a family of serine/threonine kinases that serve as targets for the small GTP binding proteins, CDC42 and RAC1, and have been implicated in a wide range of biological activities. This protein is activated by proteolytic cleavage during caspase-mediated apoptosis, and may play a role in regulating the apoptotic



events in the dying cell.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.