



Recombinant LPS-assembly protein lptD (lptD)

Product Code	CSB-BP815993EGX
Abbreviation	lptD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CWE6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	ADLASQ CMLGVPSYDR PLVQGDNDL PVTINADHAK GDYPDDAVFT GSVDIMQGNL RLQADEVQLH QKEAPGQPEP VRTVDALGNV HYDDNQVILK GPKGWANLNT KDTNVWEGDY QMVGRQGRGK ADLMKQRGEN RYTILDNGSF TSCLPGSDTW SVVGSEIHD REEQVAEIWN ARFKVGPVPI FYSPYLQLPV GDKRRSGFLI PNAKYTTTNY FEFYLPYYWN IAPNMDATIT PHYMHRRGNI MWENEFYRLS QAGAGLMELD YLPSDKVYKD EHPNDDSSRR WLFYWNHSGV MDQVWRFNVD YTKVSDPSYF NDFDNKYGSS TDGYATQKFS VGYAVQNFNA TVSTKQFQVF SEQNTSSSYSA EPQLDVNYYQ NDVGPFDTRI YGQAVHFVNT RDDMPEATRV HLEPTINLPL SNNWGSINTE AKLLATHYQQ TNLDWYNSRN TTKLAESANR VMPQFKVDGR MVFERDMEML APGYTQTLEP RAQYLYVPYR DQSKIYNYDS SLLQSDYSGL FRDRTYGGLD RIASANQVTT GVTSTRIYDDA AVERFNISVG QIYYFTESRT GDDNITWEND DKTGSLVWAG DTYWRISDRW GLRGGIQYDT RLDNVATSNS SIEYRRDEDR LVQLNYRYAS PEYIQTLPK YYSTAEQYKN GISQVGAVAS WPIADRWISIV GAYYYDTNAN KQADSMLGVQ YSSCCYAIRV GYERKLNGWD NDKQHAYVDN AIGFNIELRG LSSNYGLGTQ EMLRSNILPY QNTL
Source	Baculovirus
Target Names	lptD
Protein Names	Recommended name: LPS-assembly protein lptD Alternative name(s): Organic solvent tolerance protein
Expression Region	25-784
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.