



Recombinant Human Lymphocyte-specific protein 1 (LSP1)

Product Code	CSB-MP013214HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P33241
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEASSDPGAEEREELLGPTAQWSVEDEEEAVHEQCQHERDRQLQAQDEE GGGHVPERPKQEMLLSLKPSEAPELDEDEGFGDWSQRPEQRQQHEGAQGT LDSGEPQCRCRSEGEQEDRPGLHAYEKEDSDEVHLEELSLSKEGPGPEDTVQ DNLGAAGAEEEQEEHQKCCQPRTPSPLVLEGTIEQSSPPLSPTTKLIDRTE NRSIEKSNSVKKSQPDLPISKIDQWLEQYTQAIETAGRTPKLARQASIELPSMA VASTKSRWETGEVQAQSAAKTPSCKDIVAGDMSKKSLWEQKGGSKTSSTIKS TPSGKRYKFVATGHGKYEKVLVEGGPAP
Source	Mammalian cell
Target Names	LSP1
Protein Names	Recommended name: Lymphocyte-specific protein 1 Alternative name(s): 47 kDa actin-binding protein 52 kDa phosphoprotein Short name= pp52 Lymphocyte-specific antigen WP34
Expression Region	1-339aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length of BC001785
Target Details	This gene encodes an intracellular F-actin binding protein. The protein is expressed in lymphocytes, neutrophils, macrophages, and endothelium and may regulate neutrophil motility, adhesion to fibrinogen matrix proteins, and transendothelial migration. Alternative splicing results in multiple transcript variants encoding different isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.