



Recombinant 30S ribosomal protein S3 (rpsC)

Product Code	CSB-MP816122VFI
Abbreviation	rpsC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8DE45
Product Type	Recombinant Protein
Immunogen Species	Vibrio vulnificus (strain CMCP6)
Purity	>85% (SDS-PAGE)
Sequence	MGQKVHPNGI RLGIVKPWNA TWFANTK DFA DNLDGDFKVR QFLTKELSKA SLSRIVIERP AKSIRVTIHT ARPGVVIGKK GEDVEKLRTA VAKIAGVPAQ INIAEVRKPE LDAQLVGDSI ASQLERRVMF RRAMKRAVQN AMRLGAKGIK VEVSGRLGGA EIARSEWYRE GRVPLHTLRA DIDYATSSAH TTYGVIGIKV WIFKGEILGG MPAATEAAEP KADKPKKQRK GRK
Source	Mammalian cell
Target Names	rpsC
Protein Names	Recommended name: 30S ribosomal protein S3
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.