



# Recombinant Ribosomal RNA small subunit methyltransferase B (rsmB)

<b>Product Code</b>	CSB-MP816535EGX
<b>Abbreviation</b>	rsmB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FD12
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKQRNLRSM AAQAIEQVVE QGQSLSNILP PLQQKVSDKD KALLQELCFG VLRTLSQSDW LINKLMARPM TGKQRTVHYL IMVGLYQLLY TRIPPHAALA ETVEGAVAIK RPQLKGLING VLRQFQRQQD ELLAEFNASD ARYLHPSWLL KRLQKAYPEQ WQSIVEANNQ RPPMWLRVNR THHSRDSWLA LLDEAGMKGF PHADYPDAVQ LETPAPVHAL PGFEEGWVTV QDASAQGCMT WLAPQNGEHI LDLCAAPGGK TTHILEVAPE AQVLAVDIDE QRLSRVYDNL KRLGMKATVK QGDGRYPSQW CGEQQFDRIL LDAPCSATGV IRRHPDIKWL RRDRDIPELA QLQSEILDAI WSHLKSGGTL VYATCSVLPE ENSLQIKAFI QRTADAELCE TGTPEQPGKQ NLPGAEEGDG FFYAKLIKK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsmB
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase B EC=2.1.1.176 Alternative name(s): 16S rRNA m5C967 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RsmB
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.