



Recombinant Putrescine aminotransferase (patA)

Product Code	CSB-YP816539EGX
Abbreviation	patA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FDF5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNRLPSSASA LACSAHALNL IEKRTLDEEE MKALNREVIE YFKEHVNPGF LEYRKSVTAG GDYGAVEWQA GGLNTLVDTQ GQEFIDCLGG FGIFNVGHRN PVVVS AVQNNQ LAKQPLHSQE LLDPLRAMLA KTLAALTPGK LKYSFFCNSG TESVEAALKL AKAYQSPRGK FTFIATSGAF HGKSLGALSA TAKSTFRKPF MPLLPGFRHV PFGNIEAMLT ALNECKKTGD DVA AVILEPI QGEGGVILPP PGYLTAVRKL CDEFGALMIL DEVQTMGRT GKMFACHEHEN VQPDILCLAK ALGGGVMPIG ATIAATEEVFS VLFDNPFLHT TTFGGNPLAC AAALATINVL LEQNLPAQAE QKGDMLLDGF RQLAREYPDL VQEARGKGM LMAIEFVDNEI GYNFASEMFR QRVLVAGTLN NAKTIRIEPP LTLTIEQCEL VIKAAARKALA AMRVSVEEA
Source	Yeast
Target Names	patA
Protein Names	Recommended name: Putrescine aminotransferase EC= 2.6.1.82 Alternative name(s): Putrescine--2-oxoglutaric acid transaminase Short name= PAT Short name= PATase
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.