



# Recombinant Fatty acid oxidation complex subunit alpha (fadJ)

<b>Product Code</b>	CSB-EP816549EGX
<b>Abbreviation</b>	fadJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FFG4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEMASAFTLN VRLDNIAVIT IDVPGEKMNT LKAEFASQVR AIIKQLRENK DLRGVVFISA KPDNFIAGAD INMIANCKTA QEAEALARQG QQLMAEIHAI PVPVIAAIHG ACLGGGLELA LACHGRVCTD DPKTVLGLPE VQLGLLPGSG GTQRLPRLIG VSTALEMILT GKQLRAKQAL KLGLVDDVVP HSILLEAAVE LAKQDRPSSR PLPVRERILA GPLGRALLFK MVGKKTEQKT QGNYPATKRI LDVIETGLAQ GTSSGYDAEA RAFGELAMTP QSQUALRNIFF ASTEVKKDPG SDAPPAPLNS VGILGGGLMG GGIAYVTACK AGLPVRIKDI NPQGINHALK YSWDQLEGKV RRRHLKASER DKQLALISGT TDYRGFAHRD LIIEAVFENL ELKQQMVAEV EQNCAAHTIF ASNTSSLPIA DIAAHAARPE QVIGLHFFSP VEKMPLVEII PHAGTSAQTI ATTVKLAKKQ GKTPIVVRDK AGFYVNRILA PYINEAIRML TEGERVEHID AALVKFGFPV GPIQLLDEVG IDTGTKIIPV LEAAYGERFS APANVVSSIL NDDRKGRKNG RGFYLYGQKG RSKKKQVDPA IYPLIGAQQG GRLSAPQVAE RCVMLMLNEA VRCLDEQVIR SVRDGDIGAV FGIGFPPFLG GPFYIDSLG AGEVVAIMQR LATQYGSRFT PCERLVEMSK RGESFWKTTA TDLQ
<b>Source</b>	E.coli
<b>Target Names</b>	fadJ
<b>Protein Names</b>	Recommended name: Fatty acid oxidation complex subunit alpha Including the following 2 domains: Enoyl-CoA hydratase/3-hydroxybutyryl-CoA epimerase EC= 4.2.1.17 EC= 5.1.2.3 3-hydroxyacyl-CoA dehydrogenase EC= 1.1.1.35
<b>Expression Region</b>	1-714
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.