



Recombinant Nucleoid-associated protein YejK (yejK)

Product Code	CSB-MP816553EGX
Abbreviation	yejK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FFR9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	SLDINQIAL HQLIKRDEQN LELVLRDSSL EPTETVVEMV AELHRVYSAK NKAYGLFSEE SELAQTLRLQ RQGEEDFLAF SRAATGRLRD ELAKYPFADG GFVLFCHYRY LAVEYLLVAV LSNLSSMRVN ENLDINPTHY LDINHADIVA RIDLTEWETN PESTRYLTFL KGRVGRKQVAD FFMDFLGASE GLNAKAQNRG LLQAVDDFTA EAQLDKAERQ NVRQQVYSYC NEQLQAGEEI ELESLSKELA GVSEVSFTEF AAEKGYELEE SFPADRSTLR QLTKFAGSGG GLTINFDAML LGERIFWDPA TDTLTIKGTP PNLRDQLQRR TSGGN
Source	Mammalian cell
Target Names	yejK
Protein Names	Recommended name: Nucleoid-associated protein YejK
Expression Region	2-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.