



Recombinant Ribosomal RNA large subunit methyltransferase I (rImI)

Product Code	CSB-MP816571EGX
Abbreviation	rImI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FJ71
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSVRLVLAKG REKSLRRHP WVFSGAVARM EGKASLGETI DIVDHQGGKWL ARGAYSPASQ IRARVWTFDP SESIDIAFFT RRLQQAQKWR DWLAQKDGLD SYRLIAGESD GLPGITIDRF GNFLVLQLLS AGAEYQRAAL ISALQTLYPE CAIYDRSDVA VRKKEGMELT QGPITGELPP ALLPIEEHGM KLLVDIQHGH KTGYLDQRD SRLATRYVE NKRVLNCFYSY TGGFAVSALM GGCSQVVSVD TSHEALDIAR QNVELNKLDL SKAEFVRDDV FKLLRTYRDR GEKFDVIVMD PPKFVENKSQ LMGACRGYKD INMLAIQLLN EGGILLTFSC SGLMTSDLFQ KIIADAAIDA GRDVQFIEQF RQAADHPVIA TYPEGLYLKG FACRVM
Source	Mammalian cell
Target Names	rImI
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase I EC=2.1.1.191 Alternative name(s): 23S rRNA m5C1962 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RImI
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.