



Recombinant Endonuclease 8 (nei)

Product Code	CSB-BP816577EGX
Abbreviation	nei
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FJU5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	PEGPEIRRA ADNLEAAIKG KPLTDVWFAF PQLKSYQSRL IGQHVTHVET RGKALLTHFS NDLTLYSHNQ LYGVWRVVDV GEEPQTTRVL RVKLQTADKT ILLYSASDIE MLTPEQLTTH PFLQRVGPDV LDPNLTPEVV KERLLSPRFR NRQFAGLLLD QAFLAGLGNV LRVEILWQVG LTGNHKAKDL NAAQLDALAH ALLDIPRLSY ATRGQVDENK YHGALFRFKV FHRDGEPCER CGGIIKTTL SSRPFYWCPG CQH
Source	Baculovirus
Target Names	nei
Protein Names	Recommended name: Endonuclease 8 Alternative name(s): DNA glycosylase/AP lyase Nei EC= 3.2.2.- EC= 4.2.99.18 DNA-(apurinic or apyrimidinic site) lyase Nei Endonuclease VIII
Expression Region	2-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>