



Recombinant Carboxylate-amine ligase YbdK (ybdK)

Product Code	CSB-MP816579EGX
Abbreviation	ybdK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FK30
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MPLPDFHVSE PFTLGIELEM QVVNPPGYDL SQDSSMLIDA VKNKITAGEV KHDITESMLE LATDVCRDIN QAAGQFSAMQ KVVQLQAAADH HLEICGGGTH PFQKWQRQEV CDNERYQRTL QNFGYLIQQA TVFGQHVHVG CASGDDAIYL LHGLSRFVPH FIALSAASPY MQGTDTRFAS SRPNIFSAPP DNGPMPWVSN WQQFEALFRC LSYTTMIDSI KDLHWDIRPS PHFGTVEVRV MDTPLTLSHA VNMAGLIQAT AHWLLTERPF KHQEKDYLLY KFNRFQACRY GLEGVITDPH TGDRRSLTEA TLRLLEKITP SAHKIGASSA IEALHRQVVS GLNEAQLMRD FVADGGSLIG LVKKHCEIWA GE
Source	Mammalian cell
Target Names	ybdK
Protein Names	Recommended name: Carboxylate-amine ligase YbdK EC= 6.3.-.-
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.