



Recombinant Proline--tRNA ligase (proS)

Product Code	CSB-BP816586EGX
Abbreviation	proS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FKZ4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MRTSQYLLST LKETPADAIEV ISHQLMLRAG MIRKLASGLY TWLPTGVRVL KKVENIVREE MNNAGAIEVL MPVVQPSELW QESGRWEQYG PELLRIADRG DRPFVLGPTH EEVITDLIRN ELSSYKQLPL NFYQIQTKFR DEVRPRFGVM RSREFLMKDA YSFHTSQESL QETYDAMYAA YSKIFSRMGL DFRAVQADTG SIGGSASHEF QVLAQSGEDD VVFSDTSDYA ANIELAEAIA PKEPRAAATQ EMTLVDTPNA KTIAELVEQF NLPIEKT VKT LLVKAVEGSS FPLVALLVRG DHELNEVKAE KLPQVASPLT FATEEEIRAV VKAGPGSLGP VNMPIPVVID RTVAAMSDFA AGANIDGKHY FGINWDRDVA TPEIADIRNV VAGDPSPDGG GTLIKRGIE VGHIFQLGTK YSEALKASVQ GEDGRNQILT MGCYGIGVTR VVAAAIEQNY DERGIVWPDA IAPFQVAILP MNMHKSFRVQ ELAEKLYSEL RAQGIEVLLD DRKERPGVMF ADMELIGIPH TIVLGDRNLD NDDIEYKYRR NGEKQLIKTG DIVDYLVKQI KG
Source	Baculovirus
Target Names	proS
Protein Names	Recommended name: Proline--tRNA ligase EC= 6.1.1.15 Alternative name(s): Prolyl-tRNA synthetase Short name= ProRS
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.