



Recombinant Carnitine operon protein CaiE (caiE)

Product Code	CSB-MP816591EGX
Abbreviation	caiE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FLA7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSYYAFEGLI PVVHPTAFVH PSAVLIGDVI VGAGVYIGPL ASLRGDYGR L IVQAGANIQD GCIMHGYCDT DTIVGENGHI GHGAILHGCV IGRDALVGMN SVIMDGAVIG EESIVAAMSF VKAGFRGEKR QLLMGTPARA VRSVSDDELH WKRLNTKEYQ DLVGRCHAAL HETQPLRQME ENRPRLQGTT DVTPKR
Source	Mammalian cell
Target Names	caiE
Protein Names	Recommended name: Carnitine operon protein CaiE
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.