



Recombinant Arabidopsis thaliana Putative L-ascorbate peroxidase 6 (APX6)

Product Code	CSB-MP816763DOA
Abbreviation	APX6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GY91
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTTTASLVK TFLFRCD SFS SFKFKCKFES PAKTRLLSPA TEKHVVRSSR AWRIRCLSDD PGSSHV FVAS RRKMVLLST VQLLSHMLPQ NGNAAEIYPV MQNEIRKVVT KGKAAGVLRL VFHDAGTFEL DDHSGGINGS IAYELERPEN IGLKKSLKVL AKAKVKVDEI QPVSWADMIS VAGSEAVSIC GGPTIPVVLG RLDSAQPDPE GKLPPETLSA SGLKECFKRK GFSTQELVAL SGAHTIGSKG FGDPTVFDNA YYKILLEKPW TSTSKMTSMV GLPSDHALVQ DDECLRWVKR YAEDQDKFFE DFTNAYIKLV NSGAKWNML
Source	Mammalian cell
Target Names	APX6
Protein Names	Recommended name: Putative L-ascorbate peroxidase 6 Short name= AtAPx08 EC= 1.11.1.11
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.