



Recombinant Arabidopsis thaliana Cyclin-dependent kinase inhibitor 4 (KRP4)

Product Code	CSB-MP816765DOA
Abbreviation	KRP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GYJ3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MGKYIRKSKI DGAGAGAGGG GGGGGGGESS IALMDVVSPS SSSSLGVLTR AKSLALQQQQ QRCLLQKPSS PSSLPPTSAS PNPPSKQKMK KKQQQMND CG SYLQLRSRRL QKKPPIVVIR STKRRKQQR NETCGRNPNP RSNLDSIRGD GSRSDSVSES VVFGKDKDLI SEINKDPTFG QNFFDLEEEH TQSFNRTTRE STPCSLIRRP EIMTPGSST KLNICVSESN QREDSLRSRSH RRRPTTPEMD EFFSGAEEEQ QKQFIEKYNF DPVNEQPLPG RFEWTKVDD
Source	Mammalian cell
Target Names	KRP4
Protein Names	Recommended name: Cyclin-dependent kinase inhibitor 4 Alternative name(s): Inhibitor/interactor of CDK protein 7 KIP-related protein 4
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.