



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 64 (Os07g0566200, LOC_Os07g37890)

Product Code	CSB-EP816784OFG-B
Abbreviation	Os07g0566200, LOC_Os07g37890
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H4S6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGNCVARS GT AVDAGGDGGE DGKRRRRRWK APREDQLGMV PGRIFSNDGR SRTATVYTQQ GRKGINQDAM LVWDGFGGED DGVLCGVFDG HGPHGHVVAR RVRDSLPLRL MSAARDSGAD MPAAAWRKAF ARAYKAMD KD LRSHPSLDCF CSGSTAVTVL KLGSDLYMAN IGDSRAVLGS REATGGGMVA VQLTVDLKPD VPSEAERIKK CRGRVFALQD EPEVPRVWLP FDDAPGLAMA RAFGDFCLKD YGVISVPEFF HWSLTEKDQF VILASDGVWD VLSNQEA VDI VSASPSRSKA AKSLVEAATR EWKTKYPTSK IDDCAVVCLY LDGKMDHERD STASLDNISI EEGSVADPNE PQEQEPTLTR NFTVRTVAGS TQEKTLAGVD ARIAGVANDQ NWSGLDGVTR VNSLVQLPRF SEERAIG
Source	E.coli
Target Names	Os07g0566200
Protein Names	Recommended name: Probable protein phosphatase 2C 64 Short name= OsPP2C64 EC= 3.1.3.16
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.