



# Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 64 (Os07g0566200, LOC\_Os07g37890)

<b>Product Code</b>	CSB-MP816784OFG
<b>Abbreviation</b>	Os07g0566200, LOC_Os07g37890
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8H4S6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGNCVARS GT AVDAGGDGGE DGKRRRRRWK APREDQLGMV PGRIFSNDGR SRTATVYTQQ GRKGINQDAM LVWDGFGGED DGVLCGVFDG HGPHGHVVAR RVRDSLPLRL MSAARDSGAD MPAAAWRKAF ARAYKAMD KD LRSHPSLDCF CSGSTAVTVL KLGSDLYMAN IGDSRAVLGS REATGGGMVA VQLTVDL KPD VPSEAERIKK CRGRVFALQD EPEVPRVWLP FDDAPGLAMA RAFGDFCLKD YGVISVPEFF HWSLTEKDQF VILASDGVWD VLSNQEA VDI VSASPSRSKA AKSLVEAATR EWKTKYPTSK IDDCAVVCLY LDGKMDHERD STASLDNISI EEGSVADPNE PQEQEPTLTR NFTVRTVAGS TQEKTLAGVD ARIAGVANDQ NWSGLDGVTR VNSLVQLPRF SEERAIG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Os07g0566200
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 64 Short name= OsPP2C64 EC= 3.1.3.16
<b>Expression Region</b>	1-427
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.