



Recombinant Cebus apella Superoxide dismutase [Mn], mitochondrial (SOD2)

Product Code	CSB-BP816815DQU
Abbreviation	SOD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8HXP1
Product Type	Recombinant Protein
Immunogen Species	Sapajus apella (Brown-capped capuchin) (Cebus apella)
Purity	>85% (SDS-PAGE)
Sequence	KHSLPDLPYD YGALEPHINA QIMQLHHSKH HAAYVNNLND TEEKYKEALA KGDVTAQIAL QPALKFNGGG HINHSIFWTN LSPNGGGEPK GELLEAIKRD FGSFDKFKER LTAASVGVQG SGWGWLGFNK ERGHLQIAAC PNQDPLQGTT GLIPLLGIDV WEHAYYLQYK NVRPDYKAI WNVINWENV T ERYMACKK
Source	Baculovirus
Target Names	SOD2
Protein Names	Recommended name: Superoxide dismutase [Mn], mitochondrial EC= 1.15.1.1
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.