



Recombinant *Drosophila erecta* Endophilin-A (endoA)

Product Code	CSB-EP816842DLK
Abbreviation	endoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8I1I3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAFAGLKKQI NKANQYMTEK MGGAEGTKLD MDFMEMERKT DVTVELVEEL QLKTKEFLQP NPTARAKMAA VKGISKLSGQ AKSNTYPQPE GLLAECMLTY GKKLGEDNSV FAQALVEFGE ALKQMADV KY SLDDNIKQNF LEPLHHMQTK DLKEVMHHRK KLQGRRLDFD CKRRRQAKDD EIRGAEDKFG ESLQLAQVGM FNLLNDTEH VSQLVTFAEA LYDFHSQCAD VLRGLQETLQ EKRSEAESRP RNEFVPKTL DLNLDGGGGG LNEDGTPSHI SSSASPLPSP MRSPAKSMAV TPQRQQPCC QALYDFEPEN PGELAFKEND IITLLNRVDD NWFEGAVNGR TGYFPQSYVQ VQVPLPNGN
Source	E.coli
Target Names	EndoA
Protein Names	Recommended name: Endophilin-A Alternative name(s): SH3 domain-containing GRB2-like protein
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.