



Recombinant *Drosophila melanogaster* THO complex protein 7 (thoc7)

Product Code	CSB-EP816866DLU
Abbreviation	thoc7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IRJ8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSEQCQLRPH SDTLVRKLVE MNDEEIIKQR LLIDGDGTGE DRRIVVLLKQ FLKWASDSLD SNPIMYDRLM AQFAQCKLTA LKNVQTLQMI AGERDNYTQL VEHHEESIVL AKAEIESSKK ELITAKQIRK NKMEYDLLAS LIQDQPDSE TQRHIETIRR EIDDLVQKKL KMERKFQKRR NDFTELLMYTI HELEQQLDQD SSSSASSSSS DCDARSEPDL DDNGIMEVSD EDDDLNNSTP TKFDGARGEP KYHSVSTEDS KAMSVEEDTV LELSIDKDEH DVDVAVAN
Source	E.coli
Target Names	thoc7
Protein Names	Recommended name: THO complex protein 7
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.