



Recombinant Human Armadillo repeat-containing protein 8 (ARMC8)

Product Code	CSB-MP816877HU
Abbreviation	ARMC8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUR7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>ACLETPIR MSVLSEVTAS SRHYVDRLFDPDPQKVLQGV IDMKNAVIGN NKQKANLIVL GAVPRLLYLL QQETSSTELK TECAVVLGSL AMGTENNVKS LLDCHIIPAL LQGLLSPDLK FIEACLRLR TIFTSPVTPE ELLYTDATVI PHLMALLSRS RYTQEYICQI FSHCCKGPDH QTILFNHGAV QNIAHLLTSL SYKVRMQALK CFSVLAFENP QVSMTLNVNL VDGELLPQIF VKMLQRDKPI EMQLTSAKCL TYMCRAAIR TDDNCIVLKT LPCLVRMCSK ERLLEERVEG AETLAYLIEP DVELQRIASI TDHLIAMLAD YFKYPSVSA ITDIKRLDHD LKHAHELRLQA AFKLYASLGA NDEDIRKII ETENMMDRIV TGLSESSVKV RLAAVRCLHS LRSVQQLRT SFQDHAVWKP LMKVLQNAPD EILVVASSML CNLLLEFSPS KEPILESGAV ELLCGLTQSE NPALRVNGIW ALMNMAFQAE QKIKADILRS LSTEQLFRLL SDSDLNVLMK TLGLLRNLLS TRPHIDKIMS THGKQIMQAV TLILEGEHNI EVKEQTLCIL ANIADGTTAK DLIMTNDIL QKIKYYMGHS HVKLQLAAMF CISNLIWNEE EGSQERQDKL RDMGIVDILH KLSQSPDSNL CDKAKMALQQ YLA</p>
Source	Mammalian cell
Target Names	ARMC8
Protein Names	Recommended name: Armadillo repeat-containing protein 8
Expression Region	2-673
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.