



# Recombinant Human Armadillo repeat-containing protein 8 (ARMC8)

<b>Product Code</b>	CSB-YP816877HU
<b>Abbreviation</b>	ARMC8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IUR7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>ACLETPIR MSVLSEVTAS SRHYVDRLFDPDPQKVLQGV IDMKNAVIGN  NKQKANLIVL GAVPRLLYLL QQETSSTELK TECAVVLGSL AMGTENNVKS  LLDCHIIPAL LQGLLSPDLK FIEACLRLR TIFTSPVTPE ELLYTDATVI  PHLMALLSRS RYTQEYICQI FSHCCKGPDH QTILFNHGAV QNIAHLLTSL  SYKVRMQALK CFSVLAFENP QVSMTLVNVL VDGELLPQIF VKMLQRDKPI  EMQLTSAKCL TYMCRAAIR TDDNCIVLKT LPCLVRMCSK ERLLEERVEG  AETLAYLIEP DVELQRIASI TDHLIAMLAD YFKYPPSSVSA ITDIKRLDHD  LKHAHELRLQA AFKLYASLGA NDEDIRKKII ETENMMDRIV TGLSESSVKV  RLAAVRCLHS LRSVQQLRT SFQDHAVWKP LMKVLQNAPD EILVVASSML  CNLLLEFSPS KEPILESGAV ELLCGLTQSE NPALRVNGIW ALMNMAFQAE  QKIKADILRS LSTEQLFRLL SDSDLNVLMK TLGLLRNLLS TRPHIDKIMS  THGKQIMQAV TLILEGEHNI EVKEQTLCIL ANIADGTTAK DLIMTNDIL  QKIKYYMGHS HVKLQLAAMF CISNLIWNEE EGSQERQDKL RDMGIVDILH  KLSQSPDSNL CDKAKMALQQ YLA</p>
<b>Source</b>	Yeast
<b>Target Names</b>	ARMC8
<b>Protein Names</b>	Recommended name: Armadillo repeat-containing protein 8
<b>Expression Region</b>	2-673
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.