



Recombinant Human DNA dC->dU-editing enzyme APOBEC-3F (APOBEC3F)

Product Code	CSB-EP816878HU-B
Abbreviation	APOBEC3F
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUX4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKPHFRNTVE RMYRDTFSYN FYNRPILSRR NTVWLCYEVK TKGPSRPRLD AKIFRGQVYS QPEHHAEMCF LSWFCGNQLP AYKCFQITWF VSWTPCPDCV AKLAEFLAEH PNVTLTISAA RLYYYWERDY RRALCRLSQA GARVKIMDDE EFAYCWENFV YSEGQPFMPW YKFDDNYAFL HRTLKEILRN PMEAMYPHIF YFHFKNLRKA YGRNESWLCF TMEVVKHHSP VSWKRGVFRN QVDPETHCHA ERCFLSWFCD DILSPNTNYE VTWYTSWSPC PECAGEVAEF LARHSNVNLT IFTARLYYFW DTDYQEGLRS LSQEGASVEI MGYKDFKYCW ENFVYNDDEP FKPWKGLKYN FLFLDSKLQE ILE
Source	E.coli
Target Names	APOBEC3F
Protein Names	Recommended name: DNA dC->dU-editing enzyme APOBEC-3F EC= 3.5.4.- Alternative name(s): Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 3F
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster, thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.