



Recombinant Human Extracellular serine/threonine protein kinase FAM20C (FAM20C)

| | |
|--------------------------|---|
| Product Code | CSB-EP816901HU-B |
| Abbreviation | FAM20C |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8IXL6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | LHIALDLL PRLERRGARP SGEPGCSCAQ PAAEVAAPGW AQVRGRPGEP PAASSAAGDA GWPNKHTLRI LQDFSSDPSS NLSSHSLLEKL PPAAEPAERA LRGRDPGALR PHDPAHRPLL RDPGPRRSES PPGPGGDASL LARLFEHPLY RVAVPPLTEE DVLFNVNSDT RLSPKAAENP DWPHAGAEGA EFLSPGEEAV DSYPNWLKFH IGINRYELYS RHNPAIEALL HDLSSQRITS VAMKSGGTQL KLIMTFQNYG QALFKPMKQT REQETPPDFF YFSDYERHNA EIAAFHLDR LDFRRVPPVA GRMVNMTKEI RDVTRDKKLW RTFFISPANN ICFYGECSYY CSTEHALCGK PDQIEGSLAA FLPDLSLAKR KTWARNPWRRS YHKRKKAEWE VDPDYCEEVK QTPPYDSSHR ILDVMDMTIF DFLMGNMDRH HYETFEKFGN ETFIIHLDNG RGFYKYSHDE LSILVPLQQC CRIRKSTYLR LQLLAKEEYK LSELLMAESLR GDQVAPVLYQ PHLEALDRRL RVVLKAVRDC VERNGLHSV DDDLDEHRA ASAR |
| Source | E.coli |
| Target Names | FAM20C |
| Protein Names | Recommended name: Dentin matrix protein 4 Short name= DMP-4 Alternative name(s): Protein FAM20C |
| Expression Region | 23-584 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.