



Recombinant Human Kelch domain-containing protein 8A (KLHDC8A)

Product Code	CSB-MP816911HU
Abbreviation	KLHDC8A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYD2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVPNVKDFQ WKRLAPLPSR RVYCSLLETG GQVYAIGGCD DNGVPMDCFE VYSPEADQWT ALPRLPTARA GVAVTALGKR IMVIGGVGTN QLPLKVVEMY NIDEGKWKKR SMLREAAMGI SVTAKDYRVY AAGGMGLDLR PHNHLQHYDM LKDMWVSLAP MPTPRYAATS FLRGSKIYVL GGRQSKYAVN AFEVFDIETR SWTKFPNIPY KRAFSSFVTL DNHLYSLGGL RQGRLYRQPK FLRTMDVFDM EQGGWLKMER SFFLKKRRAD FVAGSLSGRV IVAGGLGNQP TVLETAEAFH PGKNKWEILP AMPTPRCACS SIVVKNCLLA VGGVNQGLSD AVEALCVSDS
Source	Mammalian cell
Target Names	KLHDC8A
Protein Names	Recommended name: Kelch domain-containing protein 8A Alternative name(s): Substitute for delta-EGFR expression 1 Short name= S-delta-E1
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.