



Recombinant Human Protein FAM9B (FAM9B)

Product Code	CSB-BP816923HU
Abbreviation	FAM9B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZU0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAWGKKHAG KDPVRDECEE RNRFTETREE DVTDEHGERE PFAETDEHTG ANTKKPEDTA EDLTAKRKRK KMDKTCSKTK NKSKHALRKK QLKRQKRDIYI HSLKLLNVLE EYITDEQKEE EEEEEEEEEEL IRIFQEQQKK WQQYRSVRRE RLKEMKLLRD QFVKALED FE DLCDRVFSDE DSELDN
Source	Baculovirus
Target Names	FAM9B
Protein Names	Recommended name: Protein FAM9B
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.