



Recombinant Rabbit Eukaryotic translation initiation factor 2 subunit 3 (EIF2S3)

Product Code	CSB-EP007528RB-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P33887
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥85% (SDS-PAGE)
Sequence	AGGEAGVTLG QPSLVEQDSH SGFKNELERD GGLLLIAGNE SKLKHILILG NKIDLVKESQ AKEQYGQILA FVQETVAXGA TVGQVLGAVG ALPEIFTELX IVNIGSLSTG GQVSAVKADL GKIVLTNPVX TEVGEEKSVE KHWRLIGWGQ IR
Source	E.coli
Target Names	EIF2S3
Protein Names	Recommended name: Eukaryotic translation initiation factor 2 subunit 3 Alternative name(s): Eukaryotic translation initiation factor 2 subunit gamma Short name= eIF-2-gamma
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.