



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At2g38420, mitochondrial (At2g38420)

Product Code	CSB-EP817152DOA
Abbreviation	At2g38420
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L6Y7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	VIK TLAQSSQLEN ISSVLYHLEV SEKFDTPESEI FRDVIAAYGF SGRIIEAIEV FFKIPNFRCV PSAYTLNALL LVLVRKRQSL ELVPEILVKA CRMGVRLEES TFGILIDALC RIGEVDCATE LVRYSMSQDSV IVDPRLYSRL LSSVCKHKDS SCFDVIGYLE DLRKTRFSPG LRDYTVVMRF LVEGGRGKEV VSVLNQMKCD RVEPDLVICYT IVLQGVIADE DYPKADKLFDELLLLGLAPD VYTYNVYING LCKQNDIEGA LKMMSSMNKL GSEPNVVTYN ILIKALVKAG DLSRAKTLWK EMETNGVNRN SHTFDIMISA YIEVDEVVCA HGLLEEAFNM NVFVKSSRIE EVISRLCEKG LMDQAVELLA HLV
Source	E.coli
Target Names	At2g38420
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At2g38420, mitochondrial
Expression Region	78-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.