



Recombinant Mesocricetus auratus Interleukin-12 subunit beta (IL12B)

Product Code	CSB-YP817228MRG
Abbreviation	IL12B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CJE6
Product Type	Recombinant Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Purity	>85% (SDS-PAGE)
Sequence	IWELEKDV YVVEVDWSPD AAGERVVLTC DTSEEDDIW TSDKNSEAVG SGKTLTIQVK EFSNAGQYTC HKGGKLSHS RLLHKKENG IWSTDILKDQ KDPKNKTLK CEAANYSGRF TCWWLTAIST DLKFNKSSS SSSDSRAVTC GAASLSAEKV TVDRKDYQKY SVACQEDITC PTAEETLPIG LVMEAQHKYK YENYSTGFFI RDIKPDPPK NLQLKPLRGS QMELSWEYPD SWSTPHSYFS LKFHVQVHRK RERKDESQFV DKTSATIRCS KGAEVRVRAQ DHYYNSSWSR WVSVPCS
Source	Yeast
Target Names	IL12B
Protein Names	Recommended name: Interleukin-12 subunit beta Short name= IL-12B Alternative name(s): Cytotoxic lymphocyte maturation factor 40 kDa subunit Short name= CLMF p40 IL-12 subunit p40
Expression Region	23-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.