



Recombinant 3-dehydroquinate synthase (aroB)

Product Code	CSB-BP817889EGX
Abbreviation	aroB
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FCV7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MERIVVTLGE RSYPIIASG LFNEPASFLP LKSGEQVMLV TNETLAPLYL DKVRGVLEQA GNVVDSVILP DGEQYKSLAV LDTVFTALLQ KPHGRDITLV ALGGGVVGD L TGFAAASYQR GVRFIQVPTT LLSQVDSSVG GKTAVNHPLG KNMIGAFYQP ASVVVDLDCL KTLPPRELAS GLAEVIKYGI ILDGAFNWL EENLDALLRL DGPAMAYCIR RCCELKAEVV AADERETGLR ALLNLGHTFG HAIEAEMGYG NWLHGEAVAA GMVMAARTSE RLGQFSSAET QRIITLLTRA GLPVNGPREM SAQAYLPHML RDKKVLAGE M RLILPLAIGK SEVRSGVSHE LVLNAIADCQ SA
Source	Baculovirus
Target Names	aroB
Protein Names	Recommended name: 3-dehydroquinate synthase EC= 4.2.3.4
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>