



Recombinant D-serine dehydratase (dsdA)

Product Code	CSB-EP817908EGX-B
Abbreviation	dsdA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FFF4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MENAKMNSLI AQYPLVEDLV ALKETTWFNP GTTSLAEGLP YVGLTEQDVQ DAHARLSRFA PYLAKAFPET AAAGGIIESE LVAIPAMQKR LEKEYQQPIA GQLLLKKDSH LPISGSIKAR GGIYEVLAHA EKLAL EAGLL TLEDDYSKLL SPEFKQFFSQ YSIAVGSTGN LGLSIGIMSA RIGFKVTVHM SADARAWKKA KLRSHGVTVV EYEQDYGVAV EEGRKAAQSD PNCFFIDDEN SRTLFLGYSV AGQRLKAQFA QQGRIVDADN PLFVYLPCGV GGGPGGVAFG LKLAFGDHSV CFFAEPHSP CMLLGVHTGL HDQISVQDIG IDNLTAADGL AVGRASGFVG RAMERLLDGF YTLSDQTMVD MLGWLAQEEG IRLEPSALAG MAGSQRVCAS VSYQQMHGFS AEQLRNATHL VWATGGGMVP EEEMNQYLAK GR
Source	E.coli
Target Names	dsdA
Protein Names	Recommended name: D-serine dehydratase EC= 4.3.1.18 Alternative name(s): D-serine deaminase Short name= DSD
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.