



Recombinant Galactokinase (galk)

Product Code	CSB-MP817935EGX
Abbreviation	galk
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FJS1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	SLKEKTQSL FANAFGYPAT HTIQAPGRVN LIGEHTDYND GFVLPCAIDY QTVISCAPRD DRKVRVMAAD YENQLDEFSL DAPIVAHENY QWANYVRGVV KHLQLRNNSF GGVDMMVISGN VPQGAGLSSS ASLEVAVGTV LQQLYHLPLD GAQIALNGQE AENQFVGCNC GIMDQLISAL GKGDHALLID CRSLGTKAVS MPKGVAVVII NSNFKRTLVG SEYNTRREQC ETGARFFQQP ALRDVTIEEF NAVAHELDPI VAKRVRHILT ENARTVEAAS ALEQGDLMKRM GELMAESHAS MRDDFEITVP QIDTLVEIVK AVIGDKGGVR MTGGGFGGCI VALFPEELVP AVQQAVAEQY EAKTGIKETF YVCKPSQGAG QC
Source	Mammalian cell
Target Names	galk
Protein Names	Recommended name: Galactokinase EC= 2.7.1.6 Alternative name(s): Galactose kinase
Expression Region	2-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>