



Recombinant Bifunctional protein FoID (foID)

Product Code	CSB-MP817937EGX
Abbreviation	foID
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FK41
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	AAKIIDGKT IAQQVRSEVA QKVQARIAAG LRAPGLAVVL VGSNPASQIY VASKRKACEE VGFVRSYDL PETTSEAELL ELIDALNADN TIDGILVQLP LPAGIDNVKV LERIHPDKDV DGFHPYNVGR LCQRAPRLRP CTPRGIVTLL ERYNIDTFGL NAVVIGASNI VGRPMSMELL LAGCTTTVTH RFTKNLRHHV ENADLLIVAV GKPGFIPGDW IKEGAIVIDV GINRLENGKV VGDVVFEDAA KRASYITPVP GGVGPMTVAT LIENTLQACV EYHDPQGE
Source	Mammalian cell
Target Names	foID
Protein Names	Recommended name: Bifunctional protein FoID Including the following 2 domains: Methylenetetrahydrofolate dehydrogenase EC= 1.5.1.5 Methenyltetrahydrofolate cyclohydrolase EC= 3.5.4.9
Expression Region	2-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>