



Recombinant Acetyl esterase (aes)

Product Code	CSB-YP817939EGX
Abbreviation	aes
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FK82
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKPENKLPVL DLISAEMKTV VNTLQPDLPWPATGTIAEQ RQYYTLERRF WNVGAPEMAT RAYRVPTKYG QVKTRIFYPQ PDRPATLFYL HGGGFILGNL DTHDRIMRLL ASYSQCTVIG IDYTLSPEAR FPQAIEEIVA ACCYFHQQAE DYQINMSRIG FAGDSAGAML ALASALWLRD KQIDCGKVAG VLLWYGLYGL RDSVTRLLG GVWDGLTQQD LQMYEEAYLS NDADRESPYY CLFNNDLTRE VPPCFIAGAE FDPLLDDSRLLYQTLAAHQQ PCEFKLYPGT LHAFLHYSRM MKTADALRD GAQFFTAQL
Source	Yeast
Target Names	aes
Protein Names	Recommended name: Acetyl esterase EC= 3.1.1.-
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.