



Recombinant Arabidopsis thaliana KH domain-containing protein At1g09660 (At1g09660)

Product Code	CSB-MP818113DOA
Abbreviation	At1g09660
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GWR3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MMESGAGFVA MEERISPGSF FQYPLSGFRA SPNRSPCPPS DRERYLTEL QERQKLGPF QVMPNCCRLL NHEIRRVSSF PDLDRYEHGS PFRSLGQPTN GKLDLEGWSM MQAEENCHLQ RASPFRGPSP VGWIGMPGLP NPPIVKKVIR LDVPVDKYPS YNFVGRILGP RGNSLKRVEL ATHCRVFIRG RGSVKDTVKE EKLKKGKPGYE HLCEPLHVL I EAELPEDIIN SRLEHAVHFL ESSLKPMDES MDHYKREQLK ELAALNGTLR EESPSPSLSP CLSPSMSPFN SKRAKTEI
Source	Mammalian cell
Target Names	At1g09660/At1g09670
Protein Names	Recommended name: KH domain-containing protein At1g09660 Alternative name(s): Quaking-like protein 5
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.