



Recombinant Arabidopsis thaliana GATA transcription factor 28 (GATA28)

Product Code	CSB-EP818139DOA-B
Abbreviation	GATA28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H1G0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDDLHGSNAR MHIREAQDPM HVQFEHHLH HIHNGSGMVD DQADDGNAGG MSEGVEDIP SHPGNVTDNR GEVVDRGSEQ GDQLTSLFQG QVYVFDVLP EKQAVLLLL GGRELQAAP PGLGSPHQNN RVSSLPGTPQ RFSIPQLAS LVRFREKRKG RNFDDKIRYT VRKEVALRMQ RNKGQFTSAK SNNDEAASAG SSWGNSQTWA IESSEAQHQE ISCRHCGIGE KSTPMMRRGP AGPRTLACNAC GLMWANKGAF RDLSKASPQT AQNPLPNKNE DANLETDHQI MITVANDISN SQ
Source	E.coli
Target Names	GATA28
Protein Names	Recommended name: GATA transcription factor 28 Alternative name(s): Protein TIFY 2A ZIM-like 2 protein
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.