



Recombinant *Oryza sativa* subsp. *japonica* Actin-depolymerizing factor 9 (ADF9)

Product Code	CSB-BP818142OFG
Abbreviation	ADF9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H2P8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MANSASGLAV NDECKFKFQE LKTRRGFRFI VFKIDDKAME IKVERLGQTA EGYEDFAATL PADECRYAVY DLDFVTDENC QKSKIFFFSW SPDTARTRSK MLYASSKDRF RRELDGIQCE IQATDPSEMS LDIIRARAH
Source	Baculovirus
Target Names	ADF9
Protein Names	Recommended name: Actin-depolymerizing factor 9 Short name= ADF-9 Short name= OsADF9
Expression Region	1-139
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.