



# Recombinant Human E3 ubiquitin-protein ligase SIAH1 (SIAH1)

<b>Product Code</b>	CSB-YP818230HU
<b>Abbreviation</b>	SIAH1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IUQ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTGKATPPSLYSWRGVLFTCLPAARTRKRKEMSRQTATALPTGTSKCPPSQR VPALTGTTASNNDLASLFECPCFDYVLPPIQCQSGHLVCSNCRPKLTCCPT CRGPLGSIRNLAMEKVANSVLFPCKEYASSGCEITLPHEKADHEELCEFRPYS CPCPGASCKWQGSGLDAVMPHLMHQHKSITTLQGEDIVFLATDINLPGAVDWV MMQSCFGFHFMLVLEKQEKYDGHQQFFAIVQLIGTRKQAENFAYRLELNGHR RRLTWEATPRSIHEGIATAIMNSDCLVFDTSIAQLFAENGNLGINVTISM
<b>Source</b>	Yeast
<b>Target Names</b>	SIAH1
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase SIAH1 EC= 6.3.2.-Alternative name(s): Seven in absentia homolog 1 Short name= Siah-1 Siah-1a
<b>Expression Region</b>	1-313
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length of Isoform 2
<b>Target Details</b>	This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in the development of certain forms of Parkinson s disease, the regulation of the cellular response to hypoxia and induction of apoptosis. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.